


Neonatal screening

Specimen Gate® Essentials Neonatal Screening Software Technical Requirements



Specimen Gate® Essentials is a browser based application used for the management of neonatal screening data. The application requires only a single installation to either a server or PC. Provided the application is installed to a server, other PCs on the same network with the correct security rights can access the application using a supported internet browser.

The purpose of this document is to provide a summary of the hardware and software requirements for the installation and running of the application.

SERVER REQUIREMENTS

Hardware (minimum requirements):

Quad core processor, 2.8 GHz. Memory 8 GB.

Free disk space 100 GB

Operating system:

Windows® Server 2012 R2

SQL Database:

SQL Server 2014

Web browser: Internet Explorer® 11, Firefox® v. 46 or later, Chrome™ v. 50 or later

PC REQUIREMENTS

Hardware (minimum requirements):

Quad core processor, 2.8 GHz. Memory 4 GB.

100 GB hard disk (minimum 20 GB free space required during installation)

Operating system: Windows® 7 SP1 (32/64-bit), Windows® 10 (64-bit)

SQL Database: SQL Server 2014

Web browser: Internet Explorer® 11, Firefox® v. 46 or later, Chrome™ v. 50 or later

Not for distribution in the USA or Canada. Please contact your country representative for further details.

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com

PerkinElmer, Inc.
Wallac Oy
PO Box 10
20101 Turku, Finland
Phone: + 358 2 2678 111
Fax: + 358 2 2678 357

ISO 13485
ISO 9001
CMDCAS
ISO 14001
OHSAS 18001



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2016, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

1599-9810 TCH Specimen Gate Essentials Printed in Finland